





# GET THE RELIEF YOU DESERVE

AID Ear Clip Protocol

CRANIAL ELECTROTHERAPY

1. Insert plug from ear clips into the device.
2. Saturate each felt pad with 2-3 drops of Alpha Conducting Solution™.
3. Place ear clips on ear lobes.
4. Press the power button on.
5. Increase current intensity with the “up” arrow until you feel a comfortable tingling sensation.
6. Select the timer setting for 20 if the current intensity is 3 or higher. Select the timer setting for 60 minutes if the current intensity is 1 or 2.



## Tips

- Do not treat at uncomfortably high current levels! Too much current can cause vertigo or nausea and does not lead to better results.
- Call our office for tech support if you do not tolerate an intensity setting of 1.
- It is advised to complete a CES treatment session (20 or 60 minutes) once you start it. Do not treat for only a few minutes.
- For depression, the manufacturer recommends using CES daily for 3 or 4 weeks and then 2 or 3 times per week. For anxiety and insomnia, the manufacturer recommends using CES daily for the first 1 to 2 weeks, then 2 or 3 times per week, or as needed. Plan to treat 3-8 weeks before drawing conclusions about efficacy.
- Replace the felt pads as needed. Replace batteries with alkaline AAA batteries when the battery indicator indicates no power..
- Questions? Please call us at 800-684-9343. Your Alpha-Stim® comes with tech support.



Authorized Alpha-Stim® Sales & Service  
Phone: (800) 684-9343  
Fax: (888) 684-8414  
www.AlleviaHealth.com

GET THE RELIEF YOU  
**DESERVE**



ALLEVIA  
Health

WWW.ALLEVIAHEALTH.COM

800-684-9343

**ALPHA-STIM® ORDER FORM**

Name \_\_\_\_\_

Shipping Address \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip \_\_\_\_\_

Home Phone (\_\_\_\_) \_\_\_\_\_ Day Phone (\_\_\_\_) \_\_\_\_\_ Email \_\_\_\_\_ Date of Order \_\_\_\_\_

Order at [www.AlleviaHealth.com](http://www.AlleviaHealth.com), or FAX this form to (888) 684-8414, or mail this form to:

Allevia Health, Inc.  
Authorized Alpha-Stim® Sales & Service  
922 NW Circle Blvd., Suite 160-320  
Corvallis, OR 97330

| Item Description                                                             | Item #          | Price   | Quantity | Total |
|------------------------------------------------------------------------------|-----------------|---------|----------|-------|
| <b>Alpha-Stim® M</b> for anxiety, insomnia, depression and pain (Rx only)    | ASM (400KIT)    | \$1,195 |          |       |
| <b>Alpha-Stim® AID</b> for anxiety, insomnia and depression (Rx only)        | ASAITD (500KIT) | \$795   |          |       |
| Alpha Conducting Solution™ travel size (15 ml)                               | ACS             | \$10    |          |       |
| Alpha Conducting Solution™ refill bottle (250 ml)                            | ACSR            | \$24    |          |       |
| AS-Trode™ VII Self-Adhesive Silver Electrode Pads (Pack of 4)                | AT7             | \$24    |          |       |
| AS-Trode™ VII Self-Adhesive Silver Electrode Pads (10 Packs of 4) (SAVE 37%) | AT7 (10x)       | \$150   |          |       |
| Earclip Electrode Felt Pads (Box of 256; for earclips only)                  | EEP             | \$15    |          |       |
| Probe Electrode Felt Pads (Box of 100; for probes only)                      | PEP             | \$19    |          |       |
| M Earclip Electrodes (for M model only)                                      | 401             | \$36    |          |       |
| AID Earclip Electrodes (for AID model only)                                  | 501             | \$36    |          |       |
| M Smart Probe (single probe for M model only)                                | 402             | \$15    |          |       |
| M Lead Wire (single wire for M model only)                                   | 403             | \$10    |          |       |
| M Soft Case (for M model only)                                               | 407             | \$29    |          |       |
| AID Soft Case (for AID model only)                                           | 505             | \$26    |          |       |
| M Green Lanyard (for M model only)                                           | 404             | \$5     |          |       |
| AID White Lanyard (for AID model only)                                       | 502             | \$5     |          |       |
| DVD Manual for M ( <i>Electric Body Electric Mind II.</i> )                  | 410             | \$12    |          |       |
| DVD Manual for AID                                                           | 507             | \$12    |          |       |
| Accessories for Older Alpha-Stim® Models                                     |                 |         |          |       |
| Earclip Electrodes for Alpha-Stim® 100 or Alpha-Stim® SCS                    | CE3             | \$36    |          |       |
| Set of 2 Lead Wires for Alpha-Stim® 100                                      | 2114B           | \$18    |          |       |
| Set of 2 Lead Wires for Alpha-Stim® PPM                                      | 3225            | \$18    |          |       |
| Treatment Probe for Alpha-Stim® 100 (single probe)                           | 100P            | \$15    |          |       |
| DVD Manual for Alpha-Stim® 100 ( <i>Electric Body Electric Mind</i> )        | D-EBEM          | \$29.95 |          |       |
| LED Microcurrent Tester                                                      | LED-MT          | \$10    |          |       |
| <i>Shipping: FREE for orders &gt;\$50.</i>                                   |                 |         | Subtotal |       |
| <i>For orders &lt;\$50: USPS \$5.50, UPS 2nd Day \$30, UPS Next Day \$65</i> |                 |         | Shipping |       |
| <b>5 Year Warranty &amp; Technical Support with Alpha-Stim® devices</b>      |                 |         | TOTAL    |       |

**Method of Payment**

Check enclosed (payable to Allevia Health, Inc.)

Charge my credit card:  Visa  Mastercard  AmEx  Discover

Card Number: \_\_\_\_\_ Expiration Date \_\_\_\_\_

Cardholder's Signature \_\_\_\_\_

**Return Policy:** If you are not completely satisfied with your Alpha-Stim® unit for any reason, return it within 30 days, postage paid by you, in its original salable condition, including all accessories, original packaging, manual, treatment protocols, and DVD (if purchased with Alpha-Stim® 100) for refund of the purchase price minus 15% restocking fee and minus the cost of shipping and handling at the time of purchase. Shipping charges will not be refunded. Guarantee applies only to new Alpha-Stim® stimulators purchased directly from Allevia Health, Inc. Supplies and accessories are not returnable.

# GET THE RELIEF YOU DESERVE

WWW.ALLEVIAHEALTH.COM

1-800-684-9343

## ALPHA-STIM® ACCESSORIES, SUPPLIES & TIPS



**Tip!** This large 250-ml container provides 16 refills of the travel size container (15 ml) of Alpha Conducting Solution™ that came with your device.



Alpha Conducting Solution™ (15 ml travel size)  
Item #ACS  
Models: M, AID, 100, SCS

Alpha Conducting Solution™ (250 ml refill)  
Item #ACSR  
Models: M, AID, 100, SCS

**Tip!** To prevent leakage, discard the remaining plastic ring from the cap after opening the small (15 ml) container.  
Do not discard this bottle! It is refillable. To refill, simply push the top off with your thumb.



**Tip!** Your Alpha-Stim is powered by batteries which will need replacement. There is a battery indicator on the display that monitors the battery level. For optimal freshness, purchase batteries from a large, busy retail store. We recommend high quality lithium batteries for optimal performance.



Earclip Electrode Pads  
Item #EEP  
Models: M, AID, 100, SCS

**Tip!** Remove old felt pads after you complete a CES treatment. If there is glue residue you may clean the metal surfaces of earclips with alcohol. When the surfaces are dry, place a new felt pad on each surface. The manufacturer indicates felts are intended for single use. Some clients with their own Alpha-Stim use the felts several times before replacing them.



Probe Electrode Pads  
Item #PEP  
Models: M, 100



Smart Treatment Probe  
Item #402  
Models: M only



AS-Trodes™ Self-Adhesive Silver Electrodes  
Item #AT7  
Models: M, 100, PPM

**Tip!** You can extend the life of your AS-Trodes™ by following these simple rules:

Cleanse skin with mild soap and water to remove body oils, lotions, make-up and perspiration. Rinse thoroughly to remove all soap residue and let dry before applying electrodes.

When applying electrodes, make sure there is "total contact" (no curling, gaps, etc.) with the skin.

Remove electrodes from the skin by lifting the edge of the pad, not by pulling them off by the lead wire.

Place electrodes on their protective backing and seal the bag after each treatment to prevent the gel from drying out.

If the gel becomes less sticky, apply two drops of the Alpha Conducting Solution™ to each pad to re-hydrate its surface.

When the gel loses its ability to adhere properly to the skin, it is time to replace electrodes.

If the connection between the lead wire and electrode becomes loose, you may lightly squeeze the outlet where the lead wire plugs into the electrode.




M Lead Wire (single)  
Item #403  
Models: M

Green (M model)



M Earclip Electrodes  
Item #401  
Models: M

Green (M model)



AID Earclip Electrodes  
Item #501  
Models: AID



DVD Manual for AID  
Item #507  
Models: AID only




AID Soft Case  
Item #505  
Models: AID (also fits SCS)

**Tip!** Protect your investment! This soft-sided, portable carry case is a convenient way to store your Alpha-Stim® device and its accessories.



White AID Lanyard  
Item #502  
Models: AID



DVD Manual for M  
Item #410  
Models: M only



Green M Lanyard  
Item #404  
Models: M



M Soft Case  
Item #407  
Models: M (also fits AS-100)

## ACCESSORIES FOR OLDER ALPHA-STIM® MODELS



AS-100 2 Lead Wires  
Item #2114B  
Models: 100 only

PPM 2 Lead Wires  
Item #3225  
Models: PPM only

Earclip Electrodes  
Item #CE3  
Models: 100, SCS

**Tip!** To prevent wire breakage, do not wrap earclips or leads around your Alpha-Stim device. Store all wires loosely packed in the protective case. Always remove earclips from your ear lobes by pinching the clips.

DVD Manual for AS-100  
Item #D-EBEM  
Models: 100 only

AS-100 Treatment Probe  
Item #100P  
Models: 100 only



### Can't feel current?

**Tip!** If you can not feel any tingling sensation from earclips even when the current is turned up high, check to see that the plug is all the way in, such that it is touching the unit. There should be no visible aqua color. Also, replace the battery. Otherwise, please call us for tech support. Note, most people do not feel the current when using AS-Trodes, this is normal. We recommend our LED Microcurrent Tester to test your earclips and lead wires periodically.

LED Microcurrent Tester  
for use with all Alpha-Stim models

LED Tester  
Item #LED-MT  
Models: all models